

Our screening reference

Audiometer ST 20



- Easy, fast screening test
- Optional Bone-conduction and SISI Test
- Serial PC-interface

So small and lightweight

Audiometer MA 25



- Air conduction
- Battery operated, optional AC adapter
- Continuous, pulse, warble tones
- Exceptionally lightweight, only 1,6 kg / 3.8 lb

MAICO – Innovation, Continuity, Reliability.

MAICO Product Line:

- Screening Audiometers
- Diagnostic Audiometers
- Impedance Screeners
- Middle Ear Analyzers
- OAE Test Systems
- ABR Test Systems

Professionals around the world continue to rely on MAICO products. MAICO stands for simpler, more self-explanatory operations of measuring devices, time-saving work procedures, and ergonomically optimized user interfaces.

MAICO's main focus is on the benefit for the user. MAICO is committed to state-of-the-art technology, fast reliable service and after sale support.

This brochure contains only a small segment of the comprehensive product portfolio of MAICO. To find out more about other solutions, please contact us at:

Tel.: +49 30 / 70 71 46-50
Email: export@maico.biz

Special Solutions for Children's Audiometry



Fast, reliable hearing tests in a playful manner



www.maico.biz




MAICO Diagnostic GmbH
Salzuffer 13/14,
10587 Berlin, Germany
Tel.: +49 30 / 70 71 46-50
Fax: +49 30 / 70 71 46-99
e-mail: export@maico.biz
Web: www.maico.biz

www.maico.biz



For little pilots ...

Pilot Hearing Test



"...Where is the ball?..."

The fast and fun hearing test for children

- The ideal hearing testing device for children of 2 years and older
- Select Picture Audiometry in multiple languages available
- Manual pure-tone audiometer for older children
- Large LCD display screen shows the picture of the test word

For young race driver ...

RaceCar Tymp / Audiodymp



"Win the race!"

Play-like Impedance and Reflex Screening for children

- Colorful design gives the instrument a fun look
- Animated car race keeps the child's attention
- 1000 Hz High Frequency probe tone for small children
- Full automatic test sequence - a tympanogram takes seconds
- Stapedius reflex test with 4 test frequencies

For smart detectives...

Detective / Audiometer MA 30*



The "talking" Audiometer

Manual and automatic pure tone audiometry

- Automatic Threshold Test according to Hughson Westlake
- Screening test at 30 dB in less than 2 minutes
- Frequencies freely selectable
- Stores up to 200 audiograms

* Instant print-out the results with the optional printer

Child's play, fast measurement

ERO • SCAN Screener



Fast and objective

The best solution for fast and objective hearing tests

- Frequency specific test of otoacoustic emissions with TEOAE or DPOAE
- Easy to use with PASS / REFER and graphical test result display
- Fast automatic test with accurate, reliable results
- High noise immunity allows operation in normal clinical environment
- External mini probe operates in addition to built-in probe

MAICO means fun for kids.

