

# Versatile, Portable Audiometer Series



MA 39, MA 40, MA 41 & MA 42

# MA 39, MA 40, MA 41 & MA 42



Select the model that best suits your needs - from basic air conduction to a combination of air-bone or air-bone-speech to a complete two-channel air-bone-speech audiometer with master hearing aid capabilities.

- Quiet, reliable electronic switches
- Insert phones capability
- Built-in RS232 interface for PC and NOAH compatibility
- Bright, easy to read display screen
- Lightweight, portable soft case



MA 39

MA 40

MA 41

MA 42



# Versatile, Portable Audiometers

## MA 39

The MA 39 is a portable audiometer that performs quick, accurate pure tone air conduction testing.

## MA 40

The MA 40 provides air and bone conduction testing. This unit is NOAH compatible and includes a built-in RS232 interface to transfer data to a computer.

## MA 41

The MA 41 is a portable audiometer that provides air, bone and speech testing and has the capacity to drive two free field speakers without an additional amplifier. This unit is NOAH compatible and includes a built-in RS232 interface to transfer data to a computer.

## MA 42

The MA 42 provides air, bone and speech testing as well as monaural and binaural master hearing aid testing which makes it ideal for hearing aid fitting and dispensing. It includes separate insert phone jacks, CD/tape input, monitor speaker and monitor headset. Powerful built-in amplifiers provide ample power to connect directly to your speaker system and save you the expense of buying external power amplifiers. This unit is NOAH compatible and includes a built-in RS232 interface to transfer data to a computer.

With the master hearing aid you can select four commonly used frequency slopes to approximate hearing aid acoustics and let your patients experience the difference hearing aids can make.

## Technical Data

| Part                                | MA 39 | MA 40 | MA 41 | MA 42 |
|-------------------------------------|-------|-------|-------|-------|
| <b>Standard Accessories</b>         |       |       |       |       |
| TDH 39 headset                      | ■     | ■     | ■     | ■     |
| B-71 bone receiver and cord         |       | ■     | ■     | ■     |
| Bone headband                       |       | ■     | ■     | ■     |
| Monitor/mic headset                 |       |       | ■     | ■     |
| Audiogram pad                       | ■     | ■     | ■     | ■     |
| Operation Manual                    | ■     | ■     | ■     | ■     |
| Patient response switch             |       |       |       | ■     |
| <b>Optional Accessories - Part#</b> |       |       |       |       |
| Patient response switch - 2169      | ■     | ■     | ■     |       |
| Audiocup headset - 4695             | ■     | ■     | ■     | ■     |
| Insert phones - 4790                | ■     | ■     | ■     | ■     |
| Talkback mic - 6619                 |       |       | ■     | ■     |
| Sound room kit - 6662               | ■     | ■     | ■     | ■     |
| Free field kit - 5525               |       |       | ■     | ■     |
| Sound field kit - 6608              |       |       | ■     | ■     |
| Patch cord                          | ■     | ■     | ■     | ■     |

## CPT Codes

|       |  | MA 39 | MA 40 | MA 41 | MA 42 |
|-------|--|-------|-------|-------|-------|
| 92551 | Screening, pure tone                     | ■     | ■     | ■     | ■     |
| 92552 | Threshold, air                           | ■     | ■     | ■     | ■     |
| 92553 | Air and bone                             |       | ■     | ■     | ■     |
| 92555 | Threshold, speech                        |       |       | ■     | ■     |
| 92556 | Threshold/discrimination                 |       |       | ■     | ■     |
| 92557 | Basic comprehensive                      |       |       | ■     | ■     |
| 92590 | Hearing aid exam and selection, monaural |       |       |       | ■     |
| 92591 | Binaural                                 |       |       |       | ■     |
| 92592 | Hearing aid check, monaural              |       |       |       | ■     |
| 92593 | Binaural                                 |       |       |       | ■     |

*Specifications subject to change without notice*

If air-bone-speech is air-bone-speech...  
why pay thousands more?

## MAICO Product Line:

- Screening Audiometers
- Diagnostic Audiometers
- Middle Ear Analyzers
- OAE Test Systems
- PC Compatible Products

# MAICO – Innovation, Continuity, Reliability.

This brochure contains only a small segment of the comprehensive product portfolio of MAICO. To find out more about other solutions, please contact us at:

Tel: 888-941-4201



MAICO DIAGNOSTICS  
7625 Golden Triangle Drive  
Eden Prairie, MN 55344  
Tel: 888-941-4201  
Fax: 952-903-4100  
Web: [www.maico-diagnostics.com](http://www.maico-diagnostics.com)